



Transcription for:

What Makes the RoughRider A Different Kind of Wheelchair? (Part 2)

The RoughRider's long arching foot bumpers ensure that a rider's feet remain protected at all times when bumping into any pesky obstacles. Rubber inserts also ensure that a rider's feet remain firmly planted to avoid slipping when traveling in the RoughRider. The footrests also fold to allow for easy transfer to and from the chair. Whirlwind's footrest design also allows for easy adjustments to meet any user's needs. With the use of a 13 millimeter wrench and socket, the footrest height can be changed with speed and precision.

Edited By: Allen Atristain

Narrated By: Claudia Rivas

Whirlwind Wheelchair International, 2013

www.whirlwindwheelchair.org

© 2018 Whirlwind Wheelchair International.
All rights reserved.